

Drop Shape Analyzer DSA100HP



Drop Shape Analyzer – DSA100HP
(DSA100HP1750 configuration)



Wetting and surface tension under extreme pressures and temperatures

The Drop Shape Analyzer – DSA100HP is our high-quality solution for the precise measurement of contact angle as well as surface and interfacial tension under high pressures and temperatures, primarily for tertiary oil recovery. The DSA100HP combines the DSA100 setup for optical drop shape analysis with a measuring cell for high-pressure applications up to 1750 bar. It determines the surface tension and analyzes the wetting of oil-bearing rock by surfactant solutions under the extreme pressure and temperature conditions found in reservoirs. The results help in the extraction of oil from layers of rock in the most effective manner using enhanced oil recovery (EOR) methods such as steam or surfactant flooding and thereby increase the yield.

Tasks and applications

- Surface tension of flooding solutions and their interfacial tension with oil under reservoir conditions
- Wettability and extraction of oil from oil-bearing rock
- Pressure extraction by means of liquefied gases, e.g. carbon dioxide
- Development of surfactants for liquid gases

Measuring methods and options

- Contact angle measurement using a sessile drop
- Surface tension of a liquid in gas or interfacial tension between two liquids using a pendant drop
- Upside-down pendant drop measurement of an oil drop in water
- Measurements at pressures up to 1750 bar and temperatures between -10 and 250 °C

Versatile options for dosing and measuring

With the aid of the dosing units of the DSA100HP, the drop liquid and the surrounding phase easily reach the inside of the high pressure cell. Rock samples can be positioned from the outside to measure the wettability at different positions while maintaining the pressure. The flexible design also allows contact angle measurements in a surrounding liquid phase. Precise temperature control up to 250 °C is achieved by means of an electric heater.



Sample positioning with closed pressure chamber.



Either a gas or a liquid can be added as surrounding phase.

Reliable drop shape analysis thanks to high imaging quality

The DSA100HP features a high-resolution camera and a quality zoom lens for accurate display of the drop with optimum size. The high image quality that this achieves leads to a precisely measured contact angle or surface/interfacial tension. Combined with the intelligent image evaluation algorithm of the ADVANCE software, drop shape analysis with the instrument provides exact results.

Specifications

Camera system

Connection	USB 3.0
Performance	CF04: up to 2300 fps CF06: up to 3400 fps

Optics

Zoom	7× zoom, manual
------	-----------------

Pressure control

Maximum pressure	40 to 1750 bar (580 to 25 000 psi)
------------------	------------------------------------

Temperature control

Range	up to 250 °C
-------	--------------

Pressure chamber

Material	stainless steel, Hastelloy®, or Inconel®
----------	--

Contact angle

Range	0 to 180°
Resolution	0.01°

Interfacial and surface tension

Range	0.01 to 2000 mN/m
Resolution	0.01 mN/m



Drop Shape Analyzer – DSA100 HP

Specifications



Product group specifications		DSA100HP40	DSA100HP690	DSA100HP1750
Camera system				
Connection			USB 3.0	
Resolution			CF04: 1920 x 1200 px CF06 ¹⁾ : 640 × 480 px	
Frame rate			CF04: 2300 fps CF06: 3400 fps	
Dark noise			CF04: 7 electrons CF06: 10.5 electrons	
Dynamic range			CF04: 73 dB CF06: 56.6 dB	
Optics				
Focus			manual	
Zoom			7× zoom, manual	
View angle			±4°	
Field of view			CF04: 3.9 mm × 3.9 mm to 24.7 mm × 24.7 mm CF06 ¹⁾ : 1.7 mm × 1.3 mm to 10.8 mm × 8.1 mm	
Resolution			CF04: 3.1 to 21.7 μm CF06 ¹⁾ : 2.5 to 17.8 μm	
Illumination				
Type			high power monochromatic LED	
Wave length, dominant			470 nm	
Field of light			46 mm × 46 mm (D × H)	
Dosing system				
Dosing			manual	
Drop deposition			manual	
Piston, volume		5 mL	5 mL	11 mL
Stages			x-axis ¹⁾	
Control			manual	
Length			10 mm	
Software				
ADVANCE			contact angle interfacial and surface tension of liquids	

¹⁾ optional

General specifications		DSA100HP40	DSA100HP690	DSA100HP1750
Sample dimensions				
Maximum sample space		40 cm ³	25 cm ³	25 cm ³
Maximum measuring surface		10 mm × 0 mm		
Pressure control				
Maximum pressure		40 bar (580 psi)	690 bar (10 000 psi)	1750 bar (25 000 psi)
Temperature control				
Type		fluid	electrical	electrical
Range		20 to 200 °C	20 to 200 °C	20 to 250 °C
Flow-through thermostat		yes	-	-
Temperature measurement				
Range		-50 to 400 °C		
Resolution		0.1 °C		
Precision		0.1 °C		
Accuracy		1/3 DIN B (±0.1 °C at 0 °C to ±0.8 °C at 400 °C)		
External sensor		2 connectors (PT100) ⁵⁾		
Locations		sample stage, chamber, cuvette		
Housing and peripherals				
Compartment		test liquids protected against light		
Needle protection shield		yes		
Camera und optics housing		yes		
Control keyboard		PC keyboard for ADVANCE software operation available (KB20)		
Levelling		yes		
Environment				
Operating temperature		10 to 40 °C		
Humidity		without condensation		
Instrument dimensions				
Footprint		900 mm × 600 mm (W × D)		
Height		900 mm		
Weight (without accessories)		80 kg		
Power supply				
Voltage		88 to 264 VAC		
Power consumption		100 W		
Frequency		50 to 60 Hz		
Interfaces				
PC		USB 3.0		

⁵⁾ retrofittable

Measurement specifications	DSA100HP40	DSA100HP690	DSA100HP1750
Sessile drop/Captive bubble			
Result	contact angle		
Range ²⁾	0 to 180°		
Resolution ²⁾	0.01°		
Accuracy ³⁾	0.1°		
Models	conic section, polynom, circle, Young-Laplace, height-width		
Types ⁴⁾	advancing, receding, static, dynamic		
Pendant drop/Rising drop			
Results	interfacial and surface tension		
Range	0.01 to 2000 mN/m		
Resolution	0.01 mN/m		
Model	Young-Laplace		
Types	static, dynamic		

²⁾ software-based

³⁾ instrument-based

⁴⁾ additional accessories may be required

